**Advancis**  
**Manuka honey**

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**+VISIO CARE**  
a new medical vision

**Advancis**

Medical  
[www.advancis.co.uk](http://www.advancis.co.uk)

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### Perchè Manuka honey?

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## Case Reports

### The use of Activon® honey on a non-compliant patient with chronic venous leg ulcers

Ann Smith RGN District Nursing Sister - Stocksbridge, Sheffield

**Mrs H is aged 64 years and presented to District Nurses with bilateral leg ulcers and following Doppler examination this was found to be venous in origin. Mrs H was having both legs re-dressed 3 times a week with many dressings but was non-compliant and complained of severe pain and discomfort. The patient was returning to clinic every day (twice some days) complaining of pain, offensive odour and excessive exudate.**

There were many failed attempts at treating this lady with various antimicrobials and dressings, and advice given on a healthy diet and exercise with which the patient did not comply. Mrs H's GP and nursing staff spoke to her about the possibility of surgical intervention maybe amputation of both legs, this was a shock to the patient. Fortunately we were introduced to the Activon® range by the Advancis Medical representative. The clinical staff elected to use the Activon® range as a last attempt to resolve this patient's condition. We applied the products to the patient and after two applications the odour had decreased, the exudates were slightly less and there was no pain on removal.

The patient saw an improvement and therefore began to comply with the treatment regime. As a result of this Mrs H began gentle exercise and changed her diet, which improved her condition holistically. We have continued this treatment 3 times a week for almost 6 months and have now almost healed the ulcers. The patient has now been measured for support hosiery.

Mrs H is now fully compliant with treatment, has lost 2 stones in weight and has a much healthier life style

## Activon® honey in the treatment of a fungating tumour


Caroline Farrant, Wound Care Specialist Nurse, South Wiltshire Primary Care Trust (SWPCT)

Jason is a 39 year old gentleman who was originally diagnosed with non-Hodgkin's lymphoma in 1994 following numerous operations for a non-healing abscess on his buttock which when biopsied finally confirmed the diagnosis. In the following years he had courses of radiotherapy and trips to London for various alternate treatments which also included the policies of the Bristol Cancer Centre (dietary and relaxation). He had declined any offer of chemotherapy.

Jason's disease had progressed to lymphoedema of his left arm, a large fungating tumour on his left shoulder and multiple dry cutaneous lesions on his left lateral chest wall, left arm and around his left shoulder when he was referred to the district nurse team. He had declined referral to the local hospice and palliative care team and following admission for IV antibiotics for an episode of AIE (Acute Inflammatory Episode) of his lymphoedema declared he would not have any further antibiotics.

On examination he was found to have an excess of dry scales over his anterior chest wall, areas of erythema, multiple small lesions, grade II lymphoedema and a grade II lesion with a large area of devitalized tissue, copious exudates and odour. He was pale and cachexic but denied being anaemic, the slightest exertion caused shortness of breath and he appeared weak and tired. The tumour and lymphoedema had caused obvious distortion and subluxation of his left shoulder causing physical problems with movement and dressing. Initially his treatments consisted of a hydrofibre, a carbon pad, a silver impregnated charcoal dressing, surgipads (20cm x 40cm) six a day, and tubifast (a haemostat was available if bleeding should occur) but it was obvious that these dressings did not have the capacity to contain or manage the wound exudates or odour. On consultation Jason denied pain completely, but he was never able to tolerate any compression even crepe for his lymphoedema. He stated that the odour and exudates were to him the most distressing problems.

He then requested that honey be tried. Collier (1997) suggests that it is important to ensure that individuals needs and wishes are addressed to promote autonomy and quality of life. Activon

Tulle® (supplied by Advancis Medical) was therefore applied with an absorbent pad as a secondary dressing (Eclipse® supplied by Advancis Medical).

Jason's dressings continued to be changed daily. The dry crusts were treated with Dermol 600 and Dermol 500. Jason although asked verbally denied any drawing or stinging sensation when the Activon Tulle® was applied. Despite the location of the tumour the dressing was easy to apply and easy to remove. Although observations have recorded that honey poses problems when liquefied (Lawrence 1999). Both Jason and the district nurse team scored the odour on assessment as strong despite the use of carbon. After one week the odour score was down to moderate and within the second the odour score went from slight to no odour. For Jason it enabled him to feel part of his family again and not feel ashamed. During the second week the wound had self debrided and exudates management continued to improve using Eclipse®. From six surgipads (20cm x 40cm) per day the padding was reduced to one - the Eclipse®. This was held in place using a made to measure retention vest (available from SDH) which Jason could manage despite his physical problems. No maceration was noted and it was felt without the Eclipse® dressing changes would have increased to two/three times a day.

A final thought that can only begin to sum up what it must be truly like to live with this disease. 'Can we begin to imagine what it must feel like for a patient to see part of his or her own body rotting and to have to live with the offensive smell from it, see the reaction of visitors (including doctors and nurses) and know that it signifies a lingering death'. Doyle (1980)

## The use of Algivon® on a burn to the face, neck and scalp

Jacky Edwards - Clinical Nurse Specialist - Manchester Burns Service, University Hospital of South Manchester NHS Foundation Trust

Mrs X a 34 year old woman was injured when she had applied lice lotion to her hair and then went too near to a naked flame. This caught fire and she was left with 4% Total Body Surface Area (TBSA) burns to the right side of her face, neck and scalp.

These were predominantly deep dermal in nature with a small 0.5% TBSA that was full thickness. Surgery to this area is not normally undertaken lightly as scarring can be significant, and it is difficult to recreate anatomical structures such as the ear. Initially she was treated with Atrauman and Flamazine which is standard practice, but as the wound was deeper and there was risk of more significant scarring, it was decided to treat her conservatively with Algivon®.

The patient found the dressing comfortable and soothing as well as easy to remove. The dressing was relatively easy to apply and once it had warmed up, it became very conformable and easy to use.

The dressing was changed every 2-3 days minimising nursing time and trauma to the wound bed. The promotion of autolytic debridement prevented the need for surgery, whilst maintaining an antimicrobial environment. Staff found the dressing easy to apply and remove and no adverse effects were noted.

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### References

1. Molan, P December 2001 Honey as a topical antibacterial agent for treatment of infected wounds. Available: <http://www.worldwidewounds.com/2001/November/Molan/honey-as-topical-agent.html>
2. S Natarajan, D Williamson, J Grey, KG Harding and RA Cooper (2001) Healing of an MRSA-Colonized, hydroxyurea induced leg ulcer with Honey. *Journal of Dermatological Treatment* (2001) 12, 33-66
3. Dunford C et al (2000) The use of honey in wound management. *Nursing Standard*, 15, 11, 63-68 November 2000
4. Thomas P 1998 April 1998 Odour Absorbing Dressings: A comparative laboratory study. (On-line) Available: <http://www.worldwidewounds.com/1998/march/Odour-Absorbing-Dressings/odour-absorbing-dressings.html>
5. Cooper R et al (2011) Inhibition of biofilms through the use of Manuka honey. *Wounds UK*, 2011, Vol 7, No 1
6. Kirby P et al (2009) Do honey impregnated dressings affect glycaemic control. *The Diabetic Foot Journal*, Vol 12, No 4

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**Visiocare Srl**