

Physio-Control Products and Latex

Sensitivity to medical products containing latex is an important concern in the patient care environment. In response to your inquiry, we are providing the following information about Physio-Control products. The products described below which could come in either direct or casual contact with the patient or caregiver during normal use are not manufactured with materials that contain latex as an intentionally added component or expected impurity.

The LIFEPAK® 15 Monitor/Defibrillator series and the following accessories and consumables:

- AMBU® Res-Cue Mask® Kit and Res-Cue Key® Kit
- All electrodes for Pacing, Defibrillation and ECG
- All ECG Cables
- QUIK-COMBO® Cables and Adapters
- SIGNAGEL® Electrode Gel
- Oridion Systems Ltd. Products: FilterLine® SET products for EtCO₂ circuits and intubated applications; Smart CapnoLine® and Smart CapnoLine Plus EtCO₂ circuits for non-intubated applications
- MASIMO® SET® Rainbow® SpO₂ Patient Cable and Sensors
- NIBP Hoses and Cuffs
- External Hard Paddles

The LIFEPAK 12 Defibrillator/Monitor series and the following accessories and consumables:

- AMBU Res-Cue Mask Kit and Res-Cue Key Kit
- All electrodes for Pacing, Defibrillation and ECG
- All ECG Cables
- QUIK-COMBO Cables and Adapters
- SIGNAGEL Electrode Gel
- Oridion Systems Ltd. Products: FilterLine SET products for EtCO₂ circuits and intubated applications; Smart CapnoLine and Smart CapnoLine Plus EtCO₂ circuits for non-intubated applications
- MASIMO SET SpO₂ Patient Cable and Sensors
- NELLCOR® SpO₂ Sensors and Accessories
- NIBP Hoses and Cuffs
- External Hard Paddles

Some accessories for older products no longer manufactured or distributed by Physio-Control could contain latex compounds. These accessories include the blood pressure cuffs and pressure tubing sets used with the LIFESTAT® 100 and 200 noninvasive blood pressure monitors and the VSM® 2 vital signs monitor.

The LIFEPAK 20 Defibrillator/Monitor series and the following accessories and consumables:

- AMBU Res-Cue Mask Kit and Res-Cue Key Kit
- All Electrodes for Pacing, Defibrillation and ECG
- All ECG Cables
- QUIK-COMBO Cables and Adapters
- SIGNAGEL Electrode Gel
- MASIMO SET SpO₂ Patient Cable and Sensors
- NELLCOR SpO₂ Sensors and Accessories
- External Hard Paddles

The LIFEPAK 1000 Defibrillator series and the following accessories and consumables:

- AMBU Res-Cue Mask Kit and Res-Cue Key Kit
- All Electrodes for Defibrillation and ECG
- All ECG Cables

The LIFEPAK 500 and LIFEPAK 500 DPS Automated External Defibrillators (AEDs) and the following accessories and consumables:

- AMBU Res-Cue Mask Kit and Res-Cue Key Kit
- All Electrodes for Defibrillation

The LIFEPAK CR® Plus and LIFEPAK EXPRESS® AEDs and the following accessories and consumables:

- AMBU Res-Cue Mask Kit and Res-Cue Key Kit
- All Electrodes for Defibrillation

The LUCAS™ Chest Compression System and the following accessories:

- LUCAS Air Pressure Regulator
- All Straps
- All LUCAS Connectors

For further information, please contact the Physio-Control Technical Support Team at 800.442.1142, Option 5, Option 1.

For further information, please contact Physio-Control at 800.442.1142 (U.S.), 888.879.0977 (Canada) or visit our website at www.physio-control.com.



Physio-Control Headquarters

11811 Willows Road NE
Redmond, WA 98052
www.physio-control.com

Customer Support

P. O. Box 97006
Redmond, WA 98073
Toll Free 800 442 1142
Fax 800 426 8049

Physio-Control Canada

Medtronic of Canada Ltd
99 Hereford Street
Brampton, ON
L6Y 0R3
Tel 888 879 0977
Fax 866 430 6115



Physio-Control, Inc., 11811 Willows Road NE, Redmond, WA 98052 USA